

THE UNIVERSAL 8 LT CONTAINER

Base brackets

Ideal for all requirements - can be easily secured with double sided tape to bench tops in high or low traffic areas. The 6 pack variety is versatile and can be easily incorporated into doctors and diabetics transit cases.

The 500ST offers a safe and simple alternative for personnel required to carry clinical waste containers in their glove boxes.



Mounting

IMPORTANT - the 8 litre medical waste container is designed to be used with universal brackets - wall mounted cradle, support strap or base bracket at all times for support, and is not recommended to be free standing.



Lockable plastic wall mount with a key

The 8 litre medical waste container can be installed securely by utilising the lockable plastic wall mount.

Lockable wall mounting frame

The frame can be secured to a wall and then locked closed with a padlock.

Wall straps

These wall straps enable the container to be attached to a wall but easily removed for replacement, and are ideal for non secure areas.



DISTRIBUTORS & MANUFACTURERS



Australian Factory



China Factory

P & I Sharps Containers
Telephone (08) 8294 0951
Facsimile (08) 8294 0941
Mobile 0418 824 504
Email johncook@sharpcontainers.com.au
Web www.sharpscontainers.com.au

P & I Waste Pty Ltd
ACN 008 292 582
ABN 96 008 292 582

PO Box 245
Melrose Park 5039
South Australia

Manufacturers
Australian Factory
P & I Medical Sharps Containers are manufactured by :

Blown Plastics Products Pty Ltd
(Quality Assurance License No.4053)
for P & I Sharps

Standards Australia Quality
Assurance Services Pty Ltd
ACN 050 611 642

Aust Standard
AS4031-1992



The **safest** solution for disposal of sharps, clinical and infectious waste.

CLINICAL WASTE



P&I Sharps Containers
Telephone: (08) 8294 0951
Facsimile: (08) 8294 0941
Mobile: 0418 824 504
Email: johncook@sharpcontainers.com.au
Website: www.sharpscontainers.com.au



Trolleys


Convenient trolleys are available in large sizes for the 19/35/65 litre containers and smaller sizes for the 5 and 10 litre models.







Mounting frames

These cradles enable most container to be mounted to walls or benches etc, for safety and convenience.

CONTAINERS

	VOLUME Lt	DIMENSIONS mm	DESCRIPTION lid, opening width mm	CODE	products in CARTON
	65.0	300 x 680	Yellow Standard Large Ø175mm	5PI*650 YLL	4
	19.0	260 x 415	Yellow Standard Large Ø150mm	5PI*190 YLL	8
	19.0	260 x 440	Purple Cytotoxic Large Ø150mm	5PI*190 PLL	8
	10.0	260 x 280	Yellow Standard Non spill 25 x 45mm*	5PI*100 YNS	12
	5.0	200 x 235	Yellow Standard Non spill, 25 x 45mm*	5PI*050 YNS	20
	1.4	100 x 195	Yellow Standard Non spill Ø30mm*	5PI*014 YNS	40
	0.75	100 x 130	Yellow Standard Non spill 35 x 45mm*	5PI*750YNST	40
	0.25	6pk	Yellow Standard*	5PI*0025YNS 6PK	100

*Complies with the Standards Australia Type Test 4031 License No.9004

	VOLUME Lt	DIMENSIONS mm	DESCRIPTION lid, opening width mm	CODE	products in CARTON
	35.0	300 x 450	Yellow Standard Large Ø175mm	5PI*350 YLL	4
	19.0	260 x 440	Yellow Standard Non spill 25 x 45mm*	5PI*190 YNS	8
	10.0	260 x 280	Purple Cytotoxic Non Spill 25 x 45mm*	5PI*100 PNS	12
	8.0	w 285 d 120 h 390	Yellow Standard Non spill 25 x 45mm*	5PI*080YNS	14
	2.8	155 x 210	Yellow Standard Non spill 25 x 45mm*	5PI*028 YNS	36
	1.4	90 x 210	Yellow Standard Non spill 35 x 45mm*	5PI*014YNST	40
	0.5	90 x 135	Yellow Standard Non spill 25mm*	5PI*005 YNS	40
	0.5	90 x 130	Yellow Standard Non spill 25mm*	5PI*005YNST	40

BETTER & SAFER BY DESIGN

Safety neck

The medical waste and sharps containers are fitted with a one way safety neck to prevent accidental spillage of sharps. The 19 litre is also available with a large opening snap on lid.

Built in needle notch removal system

A conveniently designed needle notch makes the removal of needles from syringes safer.

Colour coding and labelling

Colour coding is in accordance with Universal Health Authority Regulations, and the “National Guidelines for the Management of Clinical and related Wastes” as recommended by the National Health and Medical Research Council.

Labelling is clear and carries the universally accepted Bio Hazard Symbol. All P&I Sharps Medical Containers carry batch numbers.

Disposal

With the appropriate treatment requirements of the contents, the P&I Medical Sharps Containers should be incinerated in accordance with State and Federal regulations for complete and safe disposal.

Being made from Homo polymer enables the P&I Medical Sharps Containers to be incinerated without producing toxic fumes.

Storage

The design of the P&I Sharps Medical Containers enables them to be conveniently stacked side by side, and on top of each other in the carton.

Quick release straps

When a quick release strap is attached to a wall, it can provide safe mounting for containers placed on bench tops. Available for 2.8, 5, 8, 10 and 19 litre containers.



Heavy Duty Metal Syringe Safety cases

Metal cases have been designed to suit the 1.4, 8 and 65 litre P&I sharps medical waste containers.

They can be securely mounted to the floor or walls. All safety cases are constructed using double thickness steel and have twin casino locks to ensure safer environment for both staff and patrons.

Also available:

Heavy duty metal syringe safety cases in stainless steel satin finish.



CLEVER & FLEXIBLE

ACCESSORIES (available in quantities as requested)

VOLUME Lt	DESCRIPTION lid, opening width mm	CODE
1.4	Heavy Duty Metal Syringe Safety Case	5PI*014 HD
8.0	Heavy Duty Metal Syringe Safety Case	5PI*080HD
65.0	Heavy Duty Metal Syringe Safety Case	5PI*650 HD
1.4	Mounting Frame	5PI*014 FR
2.8	Mounting Frame	5PI*028 FR
5.0	Mounting Frame	5PI*050 FR
8.0	Mounting Frame	5PI*080FR
10.0	Mounting Frame	5PI*100 FR
19.0	Mounting Frame	5PI*190 FR
8.0	Lockable Plastic Wall Mount	5PI*080 WM

VOLUME Lt	DESCRIPTION lid, opening width mm	CODE
5/10	Trolley	5PI*5/10 TR
19/35/65	Trolley	5PI*19/35/65 TR
0.25	Base Bracket	5PI*0025 BB
0.5	Base Bracket	5PI*005 BB
1.4	Base Bracket	5PI*014 BB
8.0	Base Bracket	5PI*080 BB
5.0	Wall Strap	5PI*S 050 ST
8.0	Wall Strap	5PI*S 080 ST
10.0	Wall Strap	5PI*S 100 ST
19.0	Wall strap	5PI*S 190 ST